carefix[®]

Implant Stabilizer Band | حزام تثبيت غرسة الطَّعم | Pooperační prsní pás | Implantat-stabilisierungs-bånd | Kompressionsgurt für Implantate | Ζώνη σταθεροποίησης εμφυτεύματος | Banda estabilizadora de implantes | Implantin stabilointinauh | Bande stabilisante d'implant | Traka za stabiliziranje implantata | Implantátumrögzítő pánt | Fascia stabilizzatrice per impianti | インプラントスタビライザーバンド | 임플란트스테빌라이저 밴드 | Implantaatstabilisatieband | Implantatstabilisatorbånd | Passtabilizujący do implantów | Banda estabilizadora de implantes | Пояс для стабилизации имплантатов | Pás stabilizátora implantátu | Trak za stabilizaciju vsadka | Traka za stabilizaciju implantata | Bröstimplantatstöd | 植入物稳定带

Instructions for Use | تعليمات الاستعمال | Návod k použití | Brugsanvisning | Gebrauchsanweisung | Οδηγίες χρήσης | Instrucciones de uso | Käyttöohje | Instructions d'utilisation | Upute za uporabu | Használati utasítás | Istruzioni per l'uso | 使用説明 | 사용지침 | Gebruiksaanwijzing | Bruksanvisning | Instrukcja użytkowania | Instruções de uso | Инструкция по использованию | Návod na použitie | Navodila za uporabo | Uputstvo za upotrebu | Bruksanvisning | 使用说明

| | |
|----------|-------|
| EN | 4 |
| AR | 5 |
| CS | 6 |
| DA | 7 |
| DE | 8 |
| EL | 9 |
| ES | 10 |
| FI | 11 |
| FR | 12 |
| HR | 13 |
| HU | 14 |
| IT | 15 |
| JA | 16 |
| KO | 17 |
| NL | 18 |
| NO | 19 |
| PL | 20 |
| PT | 21 |
| RU | 22 |
| SK | 23 |
| SL | 24 |
| SR | 25 |
| SV | 26 |
| CN | 27 |
| ① | 28-31 |

Intended purpose

External device used to stabilize and keep implants in place after surgery.

Intended users

Women, following breast surgery.

Indications

Accessory for use with Post-Surgical Bras following breast surgery, typically augmentation.

Clinical benefits

- Keeps the implants in correct position during the healing phase

Performance characteristics**General performance characteristics**

- Breathable: Main material is breathable.
- Seamless: No seams in critical areas.

Specific performance characteristics

- Easy adjustable compression level with hook & loop closure in front
- One-size-fits-all
- Soft and stretchable material
- Lightweight

Care instructions

- Close hook & loop before washing
- Wash with similar colours
- Wash in a laundry bag
- Use of chlorine based products may decrease durability and function

Additional instructions

May be disposed with regular household waste

Only use the device after consultation and as instructed by your surgeon.

Reporting obligation

If you have any complaints, please contact your supplier directly. In the events of a serious incident that could lead to a significant deterioration of health or to death, contact the manufacturer and the competent authority of the member state. Serious incidents are defined in article 2 No. 65 of Regulation (EU) 2017/745 (MDR).

Lifetime of the device in use

The lifetime of the medical device is determined by natural wear and tear, even when treated correctly and as recommended. In case of signs of deterioration affecting performance and comfort, the device should be removed.

Contraindications

Relative contraindications:

- Very sensitive wound and skin tissue in the operated area
- Skin diseases
- Radiodermatitis
- Contact dermatitis

Warnings and precautions

- Follow the instructions for Use
- Follow your surgeon's recommendations regarding use and fitting
- Use the device according to medical indication to avoid injury or damage
- If any unexpected discomfort or wound irritation occurs, remove the device and contact your surgeon
- Make sure that the device is placed correctly to avoid pressure points and/or prevent vascular flow
- Remove the device or decrease the level of compression in case of discomfort.

Material composition

Raw materials and components used in the Post-Op Bras are tested for harmful substances according to STANDARD100 by OEKO-TEX®, product class II.

The Post-Op Bras and packaging are:

- Biologically safe, non-toxic; chlorine and dioxin free, and do not contain human or animal tissue or extracts, prohibited azo-dye pigments, medical substances or phthalates.
- Free from and unexposed to natural rubber latex or natural rubber derivatives
- Free from heavy metals or other harmful contamination and do not pose any harmful human impact associated with its intended end-use
- In compliance with REACH
- No chemicals added during treatment/finishing of the device. The device may contain softening and antistatic agents up to 1.5 %

MAIN FABRIC COMPOSITION

77% PA, 14% EL, 9% CO

① 1

Instructions for use

The device in primary packaging is prepared for direct clinical use.

Fitting the implant stabilizer band

- The device should be worn on top of a Carefix Post-Surgical Bra, high on the chest wall on the upper slope of the breast, placed just below the armpits.

| STEPS | INSTRUCTION |
|-------|--|
| 1 | N/A |
| 2 | Open the hook & loop fastener with both hands. |
| 3 | Place the device around your upper back with the hook & loop fastener corresponding to the centre of your bra. |
| 4 | N/A |
| 5 | When the device is properly placed, fasten the hook & loop above your breasts. |
| 6 | Make sure the device does not wrinkle or cause pressure in the armpit region. |

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Carefix is a registered trademark of Tytex A/S © 2022.
All rights reserved Tytex A/S, Marsvej 6, DK-7430 Ikast
Tytex A/S has been certified according to ISO 9001, ISO 14001 and ISO 13485
TYTEX A/S • Marsvej 6 • DK-7430 Ikast • www.tytex.com